


EXPIRATION EXTENSIONS FOR ANTIGEN TESTING



Extensions granted for COVID-19 antigen tests.

Where can I find the expiration date?

On the back of each test there is a "use by" date.



- Look up the new expiration date by lot number:
The lot number can be found above the "use by" date.



COMMON ANTIGEN TESTS IN ESSEX COUNTY

IF YOU HAVE SUPPLIES, PLEASE CONSIDER RE-LABELLING WITH NEW EXTENSIONS BEFORE HANDING THEM OUT SO PEOPLE KNOW THE NEW EXTENSION.

AccessBio (CareStart)

9 months from manufacture date;
3 months beyond the date originally
printed as an expiration date; and

**Example: LOT CP21J01 had an original
Exp as FEB 2022 is now MAY 2022.**



BinaxNOW (Abbott)

15 months from manufacture date;
9 months beyond the date originally
printed as an expiration date; and

**Example: LOT 134545 had an original
Exp as 16-MAY-22 (5/16/2022) is
now 16-FEB-22 (2/16/2022).**



iHealth

9 months from manufacture date (if Exp.
date is printed on or before 9/29/2022);
3 months beyond the date originally
printed as an expiration date; and

**Example: LOT 211CO21224 has a Use
BY: 2022-06-23 is now 2022-09-23.**



QuickVue (Quidel)

24 months from the manufacture date;
The date printed on the box is accurate; and
**Example: LOT F40487 Exp 2023-12-13 with
manufacture date of 2021-12-31.**

